



## ENTERALITE PUMP INTERVAL FEEDING DIRECTIONS

---

**Orders: Patient Name :** Formula: **Formula** to run at **xxx** mls /hr (rate) every **3** hours to feed a total of **xxx** mls/day.

---

- ❑ Fill the feeding bag with formula and snap the lid shut
- ❑ Load the feeding bag tubing onto the pump and shut the pump lid
- ❑ Turn on the pump, (press **ON/OFF** key and **hold for 2 seconds**)
- ❑ Program pump: Press **RATE/DOSE** key, RATE should display on pump screen. If DOSE displays on pump screen, press **RATE/DOSE** key again to switch from DOSE to RATE  
Set RATE to **xxx** ml/hr by pressing and holding down **MORE** or **LESS** key to **increase** or **decrease** the rate to desired number.
- ❑ Change RATE to DOSE on screen by pressing RATE/DOSE key. Set DOSE to **xxx** mls every 3 hrs by pressing and hold down MORE or LESS key to increase or decrease dose to desired number. DOSE increases in intervals of 10, and lowest dose is 10, not 0.
- ❑ Set INTERVAL FEEDING to **3** hrs by pressing the **FEED INT** key which will display on screen. Press **MORE** or **LESS** key to reach desired hour interval.
- ❑ You are now ready to prime the tubing. Press and hold the **PRIME** key to fill tubing with formula until it reaches the red tip.
- ❑ Flush the patient's tube with 30ml water before attaching pump tubing to patient's tube. The red tip of tubing should fit snugly.
- ❑ **START** the pump by pressing the **RUN/PAUSE** key. You know the pump is working properly when Run is displayed with arcs flashing around it on the top right corner of screen. The pump will automatically stop and start at feeding intervals set into pump's programming.
  
- ❑ You may stop the pump at any time, by pressing the **PAUSE/RUN** key and resume the feed interval by pressing the **PAUSE/RUN** key again. (Clearing the VOLUME will cause the pump to restart the feed cycle instead of continuing the feeding.)
  
- ❑ **REPLACE FEEDING SET EVERY 24 HOURS.** The pump will hold the same programming even if the pump is turned off.

To check how much formula was fed every 24 hours press the VOL/TOTAL key twice. (once shows what was fed in the last feeding interval, twice shows daily total). Press the **CLEAR** key after **VOL/TOTAL** key is pressed twice, to start a new day's feeding at zero.