



Example patient

ENTERALITE CONTINUOUS FEED PUMP DIRECTIONS

Orders: Fibersource HN 8 cans/day. Run the pump at **100 ml/ hour 20 hrs/day** to feed total volume of **2000 ml /day** as tolerated.

DIRECTIONS **Don't shake the cans of formula (prevents bubbles)**

- Fill the feeding bag with formula: **2 cans every 5 hrs (total 8 cans over 24 hrs daily)**, and snap the lid shut
- Load the feeding bag tubing onto pump and shut the pump lid
- Turn on the pump, (press **ON/OFF** key and **hold for 2 seconds**)
- Press the **RATE/DOSE** key; RATE should show on pump screen (if DOSE shows on screen, press **RATE/DOSE** key again and it will switch from DOSE to RATE)
- Set Rate to **100** ml/hr by pressing the **MORE** key. If you go too far, Press the **LESS** key to reduce number
- Press the **RATE/DOSE** key again; the **DOSE** should show on the pump screen. Set **Dose** to **2000** ml/day by pressing and holding down the **MORE** key. Once you have set the Dose and the Amount, you **will not** have to reset this again (unless the daily dose or amount of feeding changes)
- Press **RATE/DOSE** key again to make the RATE show on the pump screen.
- Press & hold down the **PRIME** key to prime the tubing with the formula. Prime the tubing until it is about one inch from the red tip end.
- Flush your feeding tube with 60ml water and connect the feeding bag by inserting the red tip into the end of your tube. It should fit snugly.
- **Start the pump** by pressing the **RUN/PAUSE** key.
- Every day when you change the feeding bag you will need to **clear out the total volume key**: After the bag change, press the **VOL TOTAL** key and then press **CLEAR**. The screen should go back to 0. Press **START/PAUSE** to restart the pump.
DO NOT START THE PUMP UNTIL YOU HAVE RESET THIS TO "0"